




STARMed Ventukit

Hood for CPAP therapy with built-in Venturi flow driver



Critical Care ▪ Patient Interfaces ▪ **STARMed**



Quality, innovation and choice

Interact with us



www.intersurgical.com



STARMed Ventukit CPAP hood with Venturi flow driver

StarMed Ventukit

StarMed Ventukit is a complete and easy to use system simply activated by connecting to the O₂ source available in the hospital departments. StarMed Ventukit unique design consists of oxygen tubes and a Venturi flow driver integral to a CPAP hood. The kit also includes an adjustable PEEP valve for immediate therapy activation.

Read more evidence on www.intersurgical.com

Technical Features

- The integrated Venturi flow driver can generate high patient flows while consuming a reduced amount of compressed oxygen
- O₂% can be adjusted from 40% to 100% at each PEEP level
- The adjustable PEEP valve can provide patient pressures in a range of settings up to 20 cmH₂O
- The user is guided by a performance table integral to the product, for independent PEEP and O₂% settings

Redesigned patient access port
With bi-directional anti-asphyxiation valve

Integrated pressure manometer

New design sealed access ports
Closer to the patient's airway

Transparent
for ease of patient monitoring

Integral Venturi flow driver
With oxygen supply tubes

Dual flow meter
Model: STS.15/30
Dual oxygen flow meter
15 + 30 L/min with O₂ coupling

Ventukit CPAP hood with Venturi flow driver, bi-directional anti-asphyxiation valve and patient access port, adjustable PEEP valve and ear plugs

Performance table subject to changes. Always refer to product labelling.

Requirement	Setting	Flow delivered
O ₂ % approx	PEEP cmH ₂ O	Flow O ₂ L/min A + B
40%	5	9 + 0
	7.5	9 + 0
	10	10 + 0
	12.5	10 + 0
	15	11 + 0
50%	20	12 + 0
	5	8 + 9
	7.5	8 + 9
	10	9 + 9
	12.5	9 + 8
	15	10 + 8
	20	12 + 9

Requirement	Setting	Flow delivered
O ₂ % approx	PEEP cmH ₂ O	Flow O ₂ L/min A + B
60%	5	7 + 16
	7.5	8 + 16
	10	8 + 14
	12.5	9 + 13
	15	9 + 13
100%	20	10 + 14
	Feed with a minimum O ₂ flow of 40 L/min	

Code	Size	Neck circumference	Box Qty.
CV100XXL/2	Extra Extra Large	≥50cm	5
CV100XL/2	Extra Large	45-52cm	5
CV100L/2	Large	40-47cm	5
CV100M/2	Medium	34-41cm	5
CV100S/2	Small	27-34cm	5
CV100XS/2	Extra Small	17-27cm	5

www.intersurgical.com/info/starmed

Make an enquiry

Dual flow meter manufactured by FLOW METER S.p.A.

Intersurgical Ltd, Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ, UK
T: +44 (0)118 965 6300 F: +44 (0)118 965 6356 info@intersurgical.com www.intersurgical.com

Management System ISO 13485:2016 ISO 9001:2015
www.tuv.com ID 900005840

Manufacturer: Intersurgical S.p.A. Mirandola (MO) Italy

Please think before you print
Save energy and paper. If you must print this information sheet please print it double sided.

IS1.21_INT • Issue 8 06.20